



ACCESALABS

Test Profile for:
33911N - eCCF 9 Panel +
Fentanyl/Meperidine/MDEA/MDMA

Substance included in test profile: Urine

<u>Substance</u>	<u>Screen</u>	<u>Confirm</u>
6-ACETYLMORPHINE (10 ng/mL SCREEN)	10	10
AMPHETAMINES (500 ng/mL SCREEN)	500	250
BARBITURATES	300	200
BENZODIAZEPINES	300	200
COCAINE METABOLITES (150 ng/mL SCREEN)	150	100
FENTANYLS (0.5 ng/mL SCREEN)	0.5	.5
MARIJUANA METABOLITES	50	15
MDA-ANALOGUES	500	250
MEPERIDINE &/OR METAB.	100	100
METHADONE	300	200
OPIATES	300	300
OXYCODONES (100 ng/mL SCREEN)	100	100
PHENCYCLIDINE	25	25
PROPOXYPHENE	300	200
CREATININE		
OXIDIZING ADULTERANTS		
pH		